




96 Well? Done!

Eppendorf *epMotion*® 96—semi-automated electronic pipette for fast and precise 96 channel microplate processing




Try it and you will see how easy and fast your 96-well and 384-well plates get processed.

»Increase your productivity and pipet whole 96-well plates at once. Just make it faster and more precise.«

The new Eppendorf epMotion 96 is an easy to use bench top system for high precision pipetting in 96- and 384-well plates. As an affordable solution it overcomes the limits of manual multichannel pipetting. Its ergonomic design and intuitive handling makes epMotion 96 a great tool for anyone who needs fast and precise liquid handling in a 96-well format. Save time and money while getting better reproducibility of your results from cell based assays, ELISA, cell seeding, nucleic acid quantification, plate reformatting or other applications.

Simultaneous 96-well pipetting

- > Performs 12 times faster than 8 channel pipets
- > Electronic pipetting with parallel piston movement for better precision and reproducibility of results
- > Precision of < 3 % at 1 μL
- > Accuracy of < 2 % at 1 μL



Ease of use

- > 0.5 to 300 μL with one system without the need to change heads
- > Auto-detection of tip size
- > Possibility to use 50 μL or 300 μL tips to maximize accuracy
- > Compact design to fit under a laminar flow hood



Ergonomics

- > Intuitive and industry proven software concept
- > Convenient touch screen
- > Easy and rapid tip load- & reloading
- > Reduced risk of repetitive stress injury (RSI)



epMotion[®] 96: One System—0.5 to 3

Compact design

Small footprint, fits under laminar flow hood and can be easily moved between labs

Control unit Apple[®] iPod[®]

High quality touch screen for convenient system control. Accessible by other iOS 7 devices (iPhone[®], iPod[®]) via WiFi

Loading frame

For quick exchange of epT.I.P.S.[®] Motion reload tip trays

Lifting table

Supports a variety of SBS-size labware such as micro plates, reservoirs, and ThermoAdapters. 384-well plates can be pipetted by indexing

Manually lift plates up and lock in place

Easy/fast adjustments to different labware heights and filling levels



00 μL **Tip attachment lever**

Minimal force needed to attach tips to the 96 nose-cone array for ergonomic operation

Auto-detect function for tip size

Convenient verification between 50 and 300 μL tips allows for more accurate volume settings and better pipetting performance

Height adjustment



Possibility to pre-set labware height for repetitive tasks or auto-pipetting mode

Accessory plate

for handling 384-well plates by indexing with the lifting table

Intuitive Operation

- > Software based on Eppendorf Xplorer® electronic pipettes provided as free Apple app
- > Intelligent, preset applications: aspiration, dilution, sequential dispensing, multi dispensing
- > Set individual parameters for aspiration or dispensing speed to match your current liquid class
- > Save frequently used volumes for quick programming
- > Possibility to program multi-step pipetting tasks for complex applications

**Applications**

- > Replication and reformatting of microplates
- > Cell seeding and media change
- > Reagent and compound addition
- > 384-wells by 4 times 96-well pipetting
- > Cell-based assays
- > ELISA handling (plate coating & wash)
- > Biochemical assays

Ease Your Workflow



- 1 Insert the appropriate tip type for 0.5–50 μL or 5–300 μL into the loading frame
- 2 Slide loading frame into epMotion 96 and apply lever to easily attach tips to the 96 nose-cone array
- 3 Select mode, e.g. »Pipette + Mix«, and adjust plate on lifting table towards the tips
- 4 Touch »Aspirate« to take liquid from source plate, change to destination plate and touch »Dispense«



Related Consumables and Accessories



300 μ L tips as reload system (color coded for easy volume identification)




50 μ L tips as reload system (reduces plastic waste of tip boxes by 40 %)




Additional loading frame




TipHolder for temporary storage of tip tray



Eppendorf twin.tec[®] PCR plates
96- and 384-well




Thermoblock for PCR 96- and 384-well
for temperature sensitive samples



Eppendorf Cell Culture Plates 96-well



Eppendorf Assay/Reader Microplates



Eppendorf Deepwell Plates

Technical Specifications

Dimension of device	Width 9 in/22.9 cm	Depth 18 in/45.7 cm	Height 21 in/53.3 cm
Weight of device	42 lb /19 kg		
Power supply	Voltage	100-240 V \pm 10 %	
	Frequency	50-60 Hz \pm 5 %	
	Max output	100 W	
Compatible plate formats	96- and 384-well microplates, shallow and deepwell (<1.3 in/3.4 cm height)		
Controller*	Capacitive touch full color LCD (via Apple iPod device WiFi connected)		
Modes:			
Pipette	Aspiration and dispensing of liquids		
Multidispense	Dispensing of liquid in multiple equal volumes		
Pipette + Mix	Aspiration and dispensing followed by user defined mixing cycles		
Manual Pipette	Aspiration of liquid can be manually stopped prior to reaching the max set volume; used for volume determination or gentle aspiration of supernatants		
Dilute + Mix	Aspiration of a concentrate and a diluent separated by an air gap. Mix after dispensing		
Multiaspirate	Aspiration of multiple equal liquid volumes; used for pooling from multiple plates		
Programming	Store sequences of multiple modes as programmed workflows		
Pipetting technology	96 channel air displacement		
Volume range	50 μ L tips	300 μ L tips	
	0.5 to 50 μ L	5 to 300 μ L	
Pipetting performance	Volume	Systematic measurement error	Random measurement error
	1 μ L (50 μ L tips)	\pm 2%	< 3%
	10 μ L (300 μ L tips)	\pm 2%	< 3%

* Apple iPod controller must be ordered separately.

Ordering information

Description	International Ordering no.	North America Ordering no.
epMotion® 96 , semi-automated electronic pipette for parallel 96 channel microplate processing (without iPod® controller)	5069 000.012	5069000004
epT.I.P.S.® Motion pipette tips		
50 μL, Reloads, Filter, PCR clean and sterile , volume range 1–50 μ L, Reloads, 12 \times 2 trays with 96 tips	0030 014.529	0030014529
50 μL, Reloads, Filter, PCR clean , volume range 1–50 μ L, Reloads, 12 \times 2 trays with 96 tips	0030 014.430	0030014430
50 μL, Reloads, Eppendorf Quality™ , volume range 1–50 μ L, Reloads, 12 \times 2 trays with 96 tips	0030 014.421	0030014421
300 μL, Reloads, Filter, PCR clean and sterile , volume range 20–300 μ L, Reloads, 12 \times 2 trays with 96 tips	0030 014.537	0030014537
300 μL, Reloads, Filter, PCR clean , volume range 20–300 μ L, Reloads, 12 \times 2 trays with 96 tips	0030 014.472	0030014472
300 μL, Reloads, Eppendorf Quality™ , volume range 20–300 μ L, Reloads, 12 \times 2 trays with 96 tips	0030 014.464	0030014464
Thermoadapter for PCR plates, 96-well , skirted, for use with 96 \times 0.2 mL tubes or a 96-well PCR plate	5075 787.008	960002199
Thermoadapter for PCR plates, 384-well , skirted, for use with a 384-well PCR plate	5075 788.004	960002202

Please visit www.eppendorf.com/epServices for available epMotion 96 Performance Plans service products.

Your local distributor: www.eppendorf.com/contact

Eppendorf AG · 22331 Hamburg · Germany
 eppendorf@eppendorf.com · www.eppendorf.com
 Eppendorf North America, Inc. · Phone: 800-645-3050
 info@eppendorf.com · www.eppendorfna.com
 Eppendorf Canada Ltd. · Phone: 800-263-8715
 canada@eppendorf.com · www.eppendorf.ca

www.eppendorf.com/automation