

Complete Mini / Mini Wide / Maxi Electrophoresis & Blotting System

Complete Mini / Mini Wide / Maxi Electrophoresis & Blotting System

The Complete Electrophoresis Systems include both modules for Gel Electrophoresis and Electroblotting and accessories, to provide a complete Mini, Mini Wide or Maxi gel casting, running and electroblotting system. The electroblotting module includes three inter-locking cassettes and six fiber pads.

A high intensity current is generated by coiled platinum electrodes. These in conjunction with the advanced cooling features allow rapid electroblotting in as little as one hour.

Features

- >Versatile: Interchangeable modules for slab gel and Electroblotting using a single universal buffer tank
- >Hinged Cassettes
- >High intensity electrodes
- >3 blot capacity
- >Rapid set-up cooling

MV-10CBS

MV-10WCBS

MV-20CBS




Ordering Information

MV-10CBS

Complete system for Mini Vertical Electrophoresis & Blotting, comprising:

- 1 x Mini Vertical Unit, including
- 2 x 2mm thick notched glass plates
- 2 x 2mm thick plain plates with 1 mm thick bounded spacers, 1 x dummy plate
- 2 x combs (1 mm thick, 12 samples), 1 x casting base, silicone mat, cooling pack




- 1x Electroblotting module, including
- 1x Internal electroblotting module
- 3x compression cassettes for gel sizes up to 100 x 100mm, 6 x fiber pads

MV-10WCBS

Complete system for Mini Wide Vertical Electrophoresis & Blotting, comprising:

- 1 x Mini Wide Vertical Unit, including
- 2 x 4mm thick notched glass plates
- 2 x 4mm thick plain plates with 1 mm thick bounded spacers, 1 x dummy plate
- 2 x combs (1 mm thick, 24 samples), 1 x casting base, silicone mat, cooling pack




- 1x Electroblotting module, including
- 1x Internal electroblotting module
- 3x compression cassettes for gel sizes up to 200 x 100mm, 6 x fiber pads

MV-20CBS

Complete system for Maxi Vertical Electrophoresis & Blotting, comprising:

- 1 x Maxi Vertical Unit including
- 2 x 4mm thick notched glass plates
- 2 x 4mm thick plain plates with 1 mm thick bounded spacers, 1 x dummy plate
- 2 x combs (1 mm thick, 24 samples), 1 x casting base, silicone mat, cooling pack



- 1x Electroblotting module including
- 1x Internal electroblotting module
- 3x compression cassettes for gel sizes up to 200 x 200mm, 6 x fiber pads



Contact Information

www.majorsci.com

info@majorsci.com

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A

T/1-408-366-9866

F/1-408-446-1107

Innovative Life Sciences Tools