

version 1.0



# 2021-2022 BROCHURE

 **major  
science**

*Innovative Life Sciences Tools*



***Power Supply***



# Power Supply



MP-310 / MP-320 / MP-510



Stackable case



Four pairs of outlet terminals



Cat. No.	MP-100	MINI-300	MINI-500	MP-310	MP-320	MP-510
		 FC CE	 FC CE <b>Upgraded</b>	 FC CE	 FC CE	 FC CE
Display	N/A	LED		2.4" TFT LCD		
Output voltage / Inc.	50V/100V	10-300V/1V	10-500V/1V	5-300V/1V		5-500V/1V
Output current / Inc.	400mA	10-400mA/1mA	1-400mA/1mA	1-700mA/1mA	10-3000mA/10mA	1-800mA/10mA
Output power / Inc.	40W	60W	120W	Output: 150W Input: 200W	300W	
Output type	Constant voltage	Constant voltage/Current		Constant voltage/Current /Power		
Timer	N/A	999 (min) with alarm/Continuous		Constant mode: 9999 (min) with alarm/Continuous Programmable mode: 999 (min) with alarm/Continuous		
Rated voltage	100V-120V~ Only			100V-240V~		
Program	N/A	2-steps (V, mA, ⌚)		Pre setting: Up to 6-steps (V, mA, W, ⌚), 30 programmed files		



\*All technical specifications are subject to change without notice.

\*Please visit our website at [www.majorsci.com](http://www.majorsci.com) for more product selection and detailed information.



#### US Office

19959 Sea Gull Way  
Saratoga, CA 95070  
U.S.A.  
T/ +1-408-366-9866  
F/ +1-408-446-1107

#### Taiwan Office

##### Head Office:

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,  
Taoyuan City 33061, Taiwan  
T/+886-3-3762878  
F/+886-3-3761310

##### Factory:

No.19, Ln. 207, Huakang St. Bade Dist.,  
Taoyuan City 33464 Taiwan  
T/ +886-3-3623319  
F/ +886-3-3623133

#### Shanghai Office

Room 612, International business exhibition center,  
9300 Hunan Road, Pudong, Shanghai, China  
National toll-free No.:400-823-9177  
T/ +86-21-50795277  
F/ +86-21-50795277

#### India Office

D. No. 12-13-99, Satguru Apts. Extn.  
Street. No. 3, Lane No.1, Tarnaka  
Secunderabad – 500 017. India  
T/ +91-40-27001515  
T/ +91-40-27001586



[www.majorsci.com](http://www.majorsci.com)

[info@majorsci.com](mailto:info@majorsci.com)